

ASLM

AFRICAN SOCIETY FOR LABORATORY MEDICINE

2018



PREVENTING AND
CONTROLLING THE
NEXT PANDEMIC:

**The Role
of the
Laboratory**


4TH INTERNATIONAL CONFERENCE
programme

“WESTGARD RULES”

10 – 13 DECEMBER 2018


TRANSCORP HILTON
Abuja, Nigeria

“Westgard Rules”




The Evolution of QC

**STEN A WESTGARD
WESTGARD QC, INC.
MADISON, WI USA
WWW.WESTGARD.COM**



THE QUEST FOR QUALITY



**How long have we running QC?
More than half a century!
(So why is it still so hard?)**



FIRST RULE: KNOW YOUR WESTGARDS

Father knows best!

"The" Westgard

- 50 years in lab medicine
- 40+ years at the University of Wisconsin
- "Westgard Rules"
- Method Validation
- Critical-Error graphs
- OPSpecs

Son knows better?

"A" Westgard

- 25 years at Westgard QC
- Publishing
- Web
- Blog
- course portal





SECOND RULE: KNOW YOUR WESTGARD WEB

Tools, Technologies and Training for Healthcare Laboratories

Westgard QC

HOME "WESTGARD RULES" ESSAYS QC APPLICATIONS LESSONS

My Cart | Check Out | Login

Beckman Coulter AU 680

Global QC Survey

HOW DO YOU **QC?**

The Great Global QC Survey

Sigma-metric analysis of Beckman Coulter AU 680.

JAMES WESTGARD
FOUNDER
Blog
About Us
Reference Materials & Resources

CALCULATORS

- QC Tools
- QC Calculators
- Method Validation Tools
- Six Sigma Calculators
- Normalized OPSpecs Calculator

HOME

Opera

Westgard On... ScienceDire... Login - PayPal AHRQ Patien... History W Managing Q... W New Year, N... Westgard O...
Web www.westgard.org/index.pl

Course Portal

Westgard QC

You must register with the website before you can enroll in a course.
You need a separate username/password in order to login at the course portal. Your westgard.com username and password will not work here.

HOME COURSES MY TRANSCRIPT MY PROFILE LOGIN LOG-OUT

Hello Visitor

Welcome to the Westgard Online Course Portal

Courses Available

New Course Available:

Risk Analysis Guidelines
An Overview

RECENT POSTS

- Westgard in the News: Advance for the Laboratory Manager, February Issue
- Coming in March: Westgard in Alexandria, Egypt
- Q & A: The Deleting Data Dilemma
- Theranos: could it get worse?
- Setting specifications for setting specifications for...
- The Great Global QC Survey
- What's New: January 2017
- Pop Quiz: What's the Sigma on Hospital Harm in the US?
- Q & A: Establishing new control ranges for unassayed chem controls
- Q & A: Is this a 4-1s violation?

RECENT ESSAYS

- FDA hears concerns about Blood Glucose Meters
- HbA1c methods in 2010, an 11 part series
- AACC 2009: Renewed Interest in Analytical Quality, but not in QC
- EQC = Eliminated QC
- The Risk Management Process

Blog:
>400 Short articles
Q&A

Website:
>60,000 members
>3 million views
>600+ essays, lessons, QC case studies, reference, resources

Course Portal:
Training in QC, Method Validation, Risk Analysis, Quality Management



The Westgard Rules

« What's New: January 2017 | Main | Setting specifications for setting specifications for... »

February 07, 2017

The Great Global QC Survey

Posted by *Sten Westgard, MS*

HOW DO YOU QC?

Tell us how YOU run your Quality Control!

There's the QC that is taught in the textbooks... and then there's the QC that's run in the real world.

We're trying to get a handle on the actual QC practices that are done around the world. So asking you to share what you do in this global survey.

To encourage you, we're offering several prizes to randomly selected survey participants. The prizes

ABOUT

Subscribe to this blog's feed

Blog powered by **Typepad**

PAYING THE BILLS

Ad closed by **Google**


Stop seeing this ad

AdChoices


WHY DO WE STILL WORRY ABOUT QC?

- **Manufacturer mean and SD used for control limits**
- **All data within 2 SD. Too good to be true!**

Control 1 Values



Control 2 Values





POOR QC = POOR PATIENT CARE

Clinical consequences of erroneous laboratory results that went unnoticed for 10 days

Tse Ping Loh, Lennie Chua Lee, Sunil Kumar Sethi et al. *J Clin Pathol* March 2013, Vol 166, No.3 260-261

- 1 test error
- 5 tests in error
- 63 results in error

Case	Primary diagnosis	Purpose of testing	Laboratory test	Units	Erroneous results	Corrected results	Potential clinical consequence	Actual clinical consequence
1	Autoimmune thyroid disease on carbimazole	Diagnostic work-up	ATG ATPO	IU/l IU/l	>3000 >1000	404 876	Repeat testing TSH: <0.02 mIU/l Free T4: 16.6 pmol/l	None
2	Syncope	Diagnostic work-up	ATG ATPO	IU/l IU/l	69 691	<20 150	None TSH: 6.90 mIU/l Free T4: 16.6 pmol/l	None
3	Partial empty sella	Disease monitoring	IGF-1	ng/ml	1509	55	Repeat testing	Repeat testing
4	Pituitary microadenoma	Disease monitoring	GH IGF-1	µg/l ng/ml	38.5 614	2.09 130	MRI imaging for suspected GH secreting adenoma	Repeat testing
5	Autoimmune thyroid disease	Disease monitoring	ATG ATPO	IU/l IU/l	96 277	<20 13	Erroneous results not seen by physician	None
6	Vitreous haemorrhage	Diagnostic work-up	ATG ATPO	IU/l IU/l	92 37	<20 <10	None TSH: 0.86 mIU/l Free T4: 16.8 pmol/l	None
7	Hypoadrenalism	Diagnostic	ACTH	pmol/l	41.1	2.1	Misdiagnosis as primary hypoadrenalism	Adrenal CT-scan ordered
8	Congenital adrenal hyperplasia	Disease monitoring	ACTH	pmol/l	102	36.6	Misdiagnosis of poor compliance to glucocorticoids	None
9	Hypothyroidism on L-thyroxine replacement	Diagnostic work-up	ATG ATPO	IU/l IU/l	126 366	23 <10	Misdiagnosis of Hashimoto's disease and need for repeat testing TSH 0.05 mIU/l Free T4: 18.4 pmol/l	None
10	Grave's disease	Diagnostic work-up	ATG ATPO	IU/l IU/l	300 >1000	<20 49	None TSH: <0.02 mIU/l Free T4: 12.7 pmol/l	None
11	Autoimmune thyroid disease	Disease monitoring	ATPO	IU/l	>1000	191	None	None
12	Hypoglycaemia for investigation	Diagnostic	GH IGF-1 Repeat testing GH IGF-1	µg/l ng/ml µg/l ng/ml	39.5 765 6.82 783	2.16 178 0.97 180	Misdiagnosis of acromegaly	None
13	Metastatic thyroid cancer	Disease monitoring	ATG	IU/l	97	<20	None	None
14	Thyroid cancer, post-surgical removal	Disease monitoring	ATG	IU/l	>3000	28	Misdiagnosis of cancer recurrence, need for further laboratory and imaging studies	None
15	Thyroid cancer, post-surgical removal	Disease monitoring	ATG	IU/l	140	<20	Misdiagnosis of cancer recurrence, need for further laboratory and imaging studies	None

The free thyroxine and thyrotropin concentrations measured together with the thyroid auto-antibody tests are provided. ACTH, adrenocorticotropic hormone (reference interval: 0.0–10.2 pmol/l), ATG, anti-thyroglobulin antibodies (negative if <40 IU/l), ATPO, anti-thyroid peroxidase antibodies (negative if <50 IU/l), GH, growth hormone (male <3.00 µg/l; female <8.00 µg/l), IGF-1, insulin-like growth factor-1 (87–238 ng/ml), free T4, free thyroxine (10.0–23.0 pmol/l), TSH, thyrotropin, (0.45–4.50 mIU/l).




THE RIGHT QC COULD HAVE CAUGHT THE ERRORS


49 patients Affected

- **4 procedures ordered in error (including CT Scan)**
- **7 patients ordered for retesting**
- **6 misdiagnoses**

Control 1 Values



Control 2 Values





HOW COULD THIS LAB MISS THIS ERROR?

CAP certified

JCI certified 2004

Singapore Service Class award 2004

ISO 15189 certified

Triple ISO certification

- ISO 9001
- ISO 14001
- ISO 18001

Awards and Awards and Awards...

We need Detailed QC HELP!





IS WASN'T JUST ONE BAD LAB – IT WAS A BAD MANUFACTURER

For 2 YEARS, Mayo Clinic: about 5% of all IGF-1 tests were false positives.

“If the Mayo Clinic observations are generalized, a laboratory performing 1000 IGF-1 tests/month would be expected to generate around 50 false-positive results each month. Some of these can be expected to lead to follow-up appointments or further testing and, ultimately, increased financial burden and anxiety for patients.”

UVA: 8-month period in 2011, “20 abnormally high IGF-1 results in 17 patients that did not agree with clinical findings. In 17 of the 20 samples, the IGF-1 concentrations measured by a mass spectrometric method were within reference intervals. In 7 of the patients, expensive growth hormone suppression tests were done; the results were within reference intervals in 6, with the result in the seventh nondiagnostic.”

Clinical Chemistry 59:8
1187–1194 (2013)

Laboratory Management

Failure of Current Laboratory Protocols to Detect Lot-to-Lot Reagent Differences: Findings and Possible Solutions

Alicia Algeciras-Schmich,¹ David E. Bruns,² James C. Boyd,² Sandra C. Bryant,³ Kristin A. La Fortune,² and Stefan K.G. Grebe^{1*}

BACKGROUND: Maintaining consistency of results over time is a challenge in laboratory medicine. Lot-to-lot reagent changes are a major threat to consistency of results.

METHODS: For the period October 2007 through July 2012, we reviewed lot validation data for each new lot of insulin-like growth factor 1 (IGF-1) reagents (Siemens Healthcare Diagnostics) at Mayo Clinic, Rochester, MN, and the University of Virginia, Charlottesville, VA. Analyses of discarded patient samples were used for comparison of lots. For the same period, we determined the distributions of reported patient results for each lot of reagents at the 2 institutions.

RESULTS: Lot-to-lot validation studies identified no reagent lot as significantly different from the preceding lot. By contrast, significant lot-to-lot changes were seen in the means and medians of 105 668 reported patient IGF-1 results during the period. The frequency of in-

allow rapid identification of between-lot result inconsistency.

© 2013 American Association for Clinical Chemistry

Maintenance of long-term stability of analytical processes and results is a pivotal task for the clinical laboratory. This process typically includes a comparison of current and new reagent lots through paired measurements of patient samples, with predefined acceptance and rejection criteria (1). Power calculations suggest that, for most assays, this approach should detect a shift in slope or intercept of 10% with 90% likelihood, if 20–30 samples are tested, provided the analytical range is not too narrow (2, 3). Each such assessment should also be compared to previous lot-to-lot evaluations to detect long-term trends. Finally, a comparison of QC values before and after a lot change, as well as external quality assurance data, might provide further data on equivalence among lots. However, in some situations,



THE GREAT GLOBAL QC SURVEY 2017



Collected March - June 2017
More than 900 labs responded
More than 105 countries represented



COUNTRIES THAT PARTICIPATED

AE - United Arab Emirates	23	ES - Spain	4	ME - Montenegro	1	RU - Russia	1
AL - Albania	1	ET - Ethiopia	3	MK - Macedonia	4	SA - Saudi Arabia	12
AM - Armenia	2	FI - Finland	1	ML - Mali	1	SD - Sudan	2
AO - Angola	2	FJ - Fiji	1	MT - Malta	1	SE - Sweden	1
AR - Argentina	5	FR - France	5	MW - Malawi	4	SG - Singapore	7
AS - American Samoa	1	FX - France, Met	3	MX - Mexico	9	SI - Slovenia	1
AT - Austria	1	GB - United Kingdom	44	MY - Malaysia	13	SV - El Salvador	1
AU - Australia	16	GH - Ghana	3	MZ - Mozambique	1	SX - Sint Maarten	1
AW - Aruba	1	GR - Greece	3	NA - Namibia	3	TH - Thailand	5
BA - Bosnia Herzegovina	8	GY - Guyana	4	NG - Nigeria	4	TR - Turkey	7
BB - Barbados	1	HK - Hong Kong	3	NL - Netherlands	6	TW - Taiwan	2
BD - Bangladesh	1	HR - Croatia	5	NO - Norway	2	TZ - Tanzania	1
BE - Belgium	7	ID - Indonesia	2	NP - Nepal	1	UA - Ukraine	2
BH - Bahrain	2	IE - Ireland	8	NZ - New Zealand	2	US - United States	350
BI - Burundi	1	IL - Israel	1	OM - Oman	6	UY - Uruguay	2
BO - Bolivia	1	IN - India	25	PA - Panama	3	VE - Venezuela	4
BR - Brazil	6	IR - Iran	4	PE - Peru	5	VI - Virgin Islands	1
BW - Botswana	3	IS - Iceland	1	PG - Papua New Guinea	1	VN - Vietnam	7
BY - Belarus	1	IT - Italy	8	PH - Philippines	10	ZA - South Africa	8
CA - Canada	41	JM - Jamaica	5	PK - Pakistan	3	ZM - Zambia	7
CH - Switzerland	3	JO - Jordan	3	PL - Poland	3	ZW - Zimbabwe	9
CL - Chile	5	JP - Japan	1	PR - Puerto Rico	9		
CN - China	3	KE - Kenya	3	PT - Portugal	10		
CO - Colombia	2	KR - Korea, South	2	QA - Qatar	5		
CR - Costa Rica	1	KW - Kuwait	6	RO - Romania	1		
CZ - Czech Republic	1	KZ - Kazakhstan	2	RS - Serbia	49		
DE - Germany	2	LB - Lebanon	12				
DK - Denmark	2	LK - Sri Lanka	2				
DZ - Algeria	1	LT - Lithuania	1				
EC - Ecuador	1	LV - Latvia	1				
EE - Estonia	2	LY - Libya	1				
EG - Egypt	6						

Top countries

US - United States	38.5%	350
RS - Serbia	5.4%	49
GB - United Kingdom	4.8%	44
CA - Canada	4.5%	41
IN - India	2.8%	25
AE - United Arab Emirates	2.5%	23






WHAT QC PRACTICES ARE BEING IMPLEMENTED IN THE “REAL WORLD”?

Anorexic QC

Gambler QC

Blind Man QC


HOW THE WORLD **QCS**




Nearly 900 laboratories
> 105 countries



ANTIQUATED QC PRACTICES




FLATTERING, BUT POSSIBLY OVERKILL




How often is the laboratory out of control? (635 labs)

30


1 of 3 labs is out-of-control every day!





HOW OFTEN DO LABS REPORT RESULTS? (EVEN THOUGH THEY'RE OUT OF CONTROL)

How often do labs over-ride the out-of-control flag and report results anyway? (648 labs)



1 out of 10 labs are ignoring their own QC!



“WESTGARD RULES”

THE ORIGINAL


**THE FIRST INNOVATION IN QUALITY
CONTROL FOR LABORATORIES**

THE ORIGINAL “WESTGARD RULES”

**Maximize error
 detection from few
 measurements**

**Attempt to balance
 work with practicality**

**Classic laboratory
 workaround**



Westgard JO, Barry PL, Hunt MR, Groth T. A multi-rule Shewhart chart for quality control in clinical chemistry. Clin Chem 1981;27:493-501.

“WESTGARD RULES”

THE PRESENT

*Are we weary and wary
of the 2s “warning rule”?*

**OPTIMIZED MULTIRULES
FOR TODAY’S INFORMATICS**





WHEN DO WE NEED “WARNING” RULES?

In the “classic/manual” multirule, the “2s warning” was used to alert operators to start checking other rules (otherwise, don’t)


Today’s labs often have QC automated by software. The computer can check all the rules all the time – no warning necessary.

In that case, what do “Westgard Rules” look like?






MODERN MULTIRULE QC PROCEDURE (N=2)



MODERN MULTIRULE QC PROCEDURE (N=3)



Use rules suited to multiples of 3

“WESTGARD RULES”

THE FUTURE

**SIX SIGMA QUALITY INTEGRATED INTO
QUALITY ASSESSMENT AND CONTROL**



THE LATEST EVOLUTION: WESTGARD SIGMA RULES



For more details:
maslablink.com
800-232-3342

WHAT ARE THE BENEFITS OF REDESIGN OF QC ?



ChiMei Hospitals, Tainan, Taiwan

- **>85% in control costs**
- **>\$50,000 annual savings in reduced reagent and control consumption**
- **>200 hours saved in troubleshooting (240 hours down to 35 hours)**

UNDERSTANDING THE IMPACT OF LABORATORY PERFORMANCE ON OUTCOMES IN A SCREENING POPULATION FOR CARDIOVASCULAR DISEASE IN TAIWAN



¹ChiMei Medical Center, Tainan, Taiwan; ²Abbott Diagnostics, Wiesbaden, Germany; ³ChiMei Medical Center, Health Management Center, Tainan, Taiwan

ISPOR 7th Asia-Pacific Conference, 3-6 September 2016, Singapore

Introduction

Risk scores for cardiovascular disease (CVD) events based on laboratory values have been established in primary prevention programs [1]. The performance of laboratory test systems may lead to discordant treatment decisions in some cases. The goal of this study was to investigate the impact of laboratory diagnostic system performance on outcomes in a screening population for CVD in Taiwan.

Methods

Data were collected from 1,396 people (Age ≥ 40 years) enrolled for CVD screening between January and April 2015 in Tainan (Table 1).

A time-to-event microsimulation model was developed (Figure 1). Starting with screening, each individual was classified into risk categories based on observed values for LDL-, HDL-, total cholesterol, and a 10-years CVD risk score. Patients with observed values of LDL \geq 190mg/L, 70<LDL<190 and a risk score \geq 7.5%, and diabetic patients with LDL between 70 and 190 plus a risk score between 5 and 7.5% were referred for treatment. They received lipid lowering drugs thus reducing risk for a CVD event. Individuals not assigned to treatment remained in the "No treatment" state until the next screening cycle after 1-4 years, the occurrence of a CVD event, or death.

Minimum and optimum test specifications as suggested in literature were tested against a control scenario assuming perfect performances (Table 2).


Samples were bootstrapped from the cohort with 100,000 iterations. Model followed a lifetime horizon and a health system perspective. Results were expressed in costs, quality

Strategy	Total-C		HDL-C		LDL-C		Source
	Bias, %	CV, %	Bias, %	CV, %	Bias, %	CV, %	
Minimum (MIN)	6.2	4.5	8.4	5.5	8.2	5.9	[16,17]
Optimum (OPT)	2.1	1.5	2.8	1.8	2.7	2.9	
Control (CON)	0.0	0.0	0.0	0.0	0.0	0.0	Assumption

Table 1. Characteristics of screening cohort (n=1,396)

Variable	Mean (SE)
Sex (% women)	35.9
BMI	24.2 (0.1)
Age, years*	52.9 (49-85)
Pulse	72.9 (0.3)
BPSys (mmHg)	116.5 (0.5)
BPDia (mmHg)	74.2 (0.3)
GLUC_FAST	100.8 (0.8)
TCHOL (mg/dL)	202.2 (1.0)
HDL-C (mg/dL)	54.7 (0.4)
LDL-C (mg/dL)	132.1 (0.9)

Figure 1. Time-to-event microsimulation model structure



Setting, risk & events	Value	Distribution	Source
CVD risk (10 years)	Risk based equations		[1]
In-hospital mortality from stroke	10.1	Beta	[2]
In-hospital mortality from AMI	6.5	Beta	[3]
Mortality from stroke at 1 year, %	12.0	Weibull	[2]
Mortality from MI at 1 year, %	6.0	Weibull	[4]
Mortality (non-CVD)	Age-, sex-specific/lifeable	Weibull	[5]
Annual risk for recurrent CVD, %	6.8	Weibull	[4]
Risk reduction under treatment, % (95%CI)	65 (58, 73)	Uniform	[6]
CVD event type (Stroke vs. MI), %	76	Beta	[7]
Smoking prevalence Male, % (95%CI)	31 (28.0, 35.2)	Uniform	[8]
Smoking prevalence Female, % (95%CI)	3.4 (2.8, 4.2)	Uniform	[8]
Laboratory results	Mean from cohort	Normal	[8]
Screening cycle, years	1 to 4	Uniform	Assumption
Utility			
Baseline	0.936	Beta	[9]
Lipid lowering treatment, Mean (SD)	0.934 (0.001)	Beta	[10]
MI (disability), Mean (SD)	0.080 (0.048)	Beta	[9]
Stroke (disutility), Mean (SD)	0.242 (0.039)	Beta	[9]
Post MI	0.799 (0.010)	Beta	[9]
Post stroke	0.576 (0.010)	Beta	[9]
Costs, NTS²			
Screening & visit	5.585	Uniform	[11]
Permanent lipid lowering treatment	15.639	Uniform	[12]
Stroke	74.832	Uniform	[13]
Post-stroke	45.630	Uniform	[14]
MI	189.497	Uniform	[15]
Post MI	82.897	Uniform	[14]
Annual discount rate for costs and utility, %	3.0		

Table 3. Model input assumptions.

(1) All costs adjusted for inflation with a 3% rate to 2015 New Taiwan dollar. (2) A range of $\pm 25\%$ was used to create upper and lower bounds.

Results


Results in terms of incremental values compared to the control scenario are summarized in Table 2.

Analytical measurement uncertainty caused by CV and bias resulted in discordant management in some cases. The MIN and OPT strategy led to different decisions in 14.1% and 4.1%, respectively. Patients who had not received preventive treatment based on erroneous results had a higher risk for CVD events at an earlier time. The observed strategies showed a small but significant per 1,000 subjects (95%CI 43, 220), whereas a non-significant trend was observed for the OPT performance. Loss in QALY resulting from unnecessary treatment, earlier and increased risk for events, or increased mortality was found to be significant for MIN but not for OPT.

CVD related life-time costs were significantly higher compared to CON for MIN (+NT\$ 8,753) and OPT (+NT\$ 2,075).

Outcome value	MIN vs. CON		OPT vs. CON		MIN vs. OPT	
	Mean	(95%CI)	Mean	95%CI	Mean	95%CI
Δ Costs per patient, NTS	8753	(7516; 9990)	2075	(844; 3307)	6678	

Figure 2. Incremental costs and QALY per strategy compared to the Control.



Microsimulation with 100,000 samples. Mean, 95%CI of Δ Costs per patient, and Δ QALY per 1,000 subjects.

The "Minimum" and "Optimum" strategy led to higher costs compared to the Control for 32% and 27%, respectively (Fig. 3).

As revealed from a sensitivity analysis, for each increase in percent point of CV, negative or positive bias 22, 43 or 7 individuals per 1,000 screened subjects would be either over- or under-treated compared to the control (Fig. 4, Tab. 3). Negative bias particularly increased the risk for denying preventive treatment, and would affect six times more patients than positive bias.

Incremental costs per patients caused from discordant management decisions and related consequences would accrue to NT\$786 (96%CI 734,838), NT\$762 (612,912) and NT\$699 (668,729), per percent increase in negative bias, positive bias and

Variable	Response	Coeff.	(95%CI)	R-Sq, %
(-) bias	Δ UT, per 1,000 subjects	43.3	(40.6; 45.9)	
Δ Costs, NTS		786	(734,838)	

Figure 3. Distribution of incremental costs (IC) per strategy.





Figure 4. Impact of increasing Bias and CV on discordant treatments. UT: Over-treated; UT: Under-treated. No significant deviation from the Control strategy was observed for both, the number of individuals over-treated with increasing negative bias, and the number of individuals under-treated with increasing positive bias.



Conclusions

- Analytical measurement uncertainty may impose a higher risk for missing prevention opportunities.
- The selection of high performance diagnostic systems plus a strict quality control management in the laboratory conforming to the optimum specification is critical to consistently providing high and efficient quality of care.

Limitations

Results are limited to the information derived from the cohort. Risk scores may not accurately estimate the actual risk. Patient characteristics were sampled from distributions in order to reflect variability and uncertainty. The model assumed a stable lipid status over the time span of the simulation.

Acknowledgement

We thank Roger Low and Sten Westgard for motivation and fruitful discussions.

References

- Stone NJ, Robinson JG, Lichtenstein AH, et al. 2013 ACC/AHA guideline on the treatment of blood cholesterol to reduce atherosclerotic cardiovascular risk in adults. *Circulation*. 2014;129(25 Suppl 2):S1-45.
- Lee HC, Chang KC, Huang YC, et al. Readmission, mortality, and first-year medical costs after stroke. *JAMA*. 2013;76(12):703-714.
- Yu WH, Lu TH, Chen NG, et al. The temporal trends of incidence, disability, and in-hospital mortality of acute myocardial infarction over 15years in a Taiwanese population. *Int J Cardiol*. 2016;209:103-113.
- Chiang FT, Shyu KG, Wu CJ, et al. Predictors of 1-year outcomes in the Taiwan Acute Coronary Syndrome Full Spectrum Registry. *J Formos Med Assoc*. 2014;113(11):794-802.
- Ministry of Interior, Department of Statistics, Taiwan. LifeTable 2013. <http://sof.w.moi.gov.tw/stat/english/elt/elt02210.htm>. Accessed May 18, 2016.
- Yaylor F, Hu H, et al. Meta-analysis of primary prevention of cardiovascular disease. *The Cochrane database of systematic reviews*. 2013(1).
- Cheng GL, Chien HC, Lee CH, Lin SJ, Yang YH. Validity of in-hospital mortality data among patients with acute myocardial infarction or stroke in the National Health Insurance Research Database in Taiwan. *Int J Cardiol*. 2015;201:36-101.
- Ng M, Freeman MK, Fleming TD, et al. Smoking prevalence and cigarette consumption in 187 countries, 1980-2012. *JAMA*. 2014; 311(2):183-92.
- Kang EJ, Ko SO, et al. The impact of tobacco and alcohol consumption on noninstitutionalized community residents in Korea. *Value Health*. 2009;12 Suppl 3:S114-117.
- Pandya A, Sy S, Cho S, Weinstein MC, Gaziano TA. Cost-effectiveness of 10-Year Risk Thresholds for Initiation of Statin Therapy for Primary Prevention of Cardiovascular Disease. *JAMA*. 2015;314(2):142-150.
- Lin YK, Chen CP, Tsai WC, Chiao YC, Lin BY. Cost-effectiveness of clinical pathway in coronary artery bypass surgery. *J Med Syst*. 2011;35(2):203-213.
- Cobiac LJ, Mosen L, Wendt J, Carter R, Vos T, et al. Improving the cardiovascular disease prevention in Australia: a modelling study. *BMC Public Health*. 2012;12(389).
- Chen HW, Chou SY, Liu JT. Hospital ownership and performance: evidence from stroke and cardiac treatment in Taiwan. *J Health Econ*. 2008;27(5):1208-1223.
- Chung CW, Wang JD, Yu CF, Yang MC. Lifetime medical expenditure and life expectancy lost attributable to smoking through major smoking related diseases in Taiwan. *Tob Control*. 2007;16(6):394-399.



HEALTH ECONOMICS OUTCOMES OF OPTIMAL SIX SIGMA QUALITY

HEOR Focus: impact of individual risk categorization and 10-year CVD score. Samples bootstrapped from [historical database] cohort with 100,000 iterations. Model followed a lifetime horizon and a health system perspective.

Tests included:

- LDL
- HDL
- total cholesterol

Variables studied:

- Minimum (low Sigma) test performance
- Optimum (high Sigma) test performance

Outcomes assessed:

- Costs of patient care
- Over- and under-treatment of patient
- Quality-adjusted life-years (QALY) of patient



HOW SIX SIGMA METHODS IMPACT QC AND PATIENT OUTCOMES

Discordant Treatment Impact:

For each percentage increase in...

+1% Imprecision = + 22 / 1,000 patients over/under-treated

-1% Negative bias = + 43 / 1,000 patients undertreated

+1% Positive bias = + 22 / 1,000 patients overtreated

Cost Impact:


For each percentage increase in...

-1% Negative bias = **NT\$786 (96%CI 734;838),**

+1% Negative bias = **NT\$762 (612;912)**

+1% Imprecision = **NT\$699 (668;729)**

Discordant treatments per 1,000 subjects



HOW SIX SIGMA METHODS IMPACT QC AND PATIENT OUTCOMES

Impact of both CV and Bias on Discordant Management:

Minimum (low Sigma) causes 14.1% of patients to incur discordant health management

Optimum (high Sigma) causes 4.1% of patients to incur discordant health management

LOSS OF LIFE YEARS PER 1,000 PATIENTS:

Minimum (low Sigma) causes 131 Life Years Loss


Optimum (high Sigma) causes no statistically significant loss of life versus perfect scenario.


CVD LIFETIME COSTS PER PATIENT:

MIN (+NT\$ 8,753) vs. OPT (+NT\$ 2,075).

Figure 2. Incremental costs and QALY per strategy compared to the Control.

Microsimulation with 100,000 samples. Mean, 95%CI of Δ Costs per patient, and Δ QALY per 1,000 subjects.





CONCLUSION

QC must evolve with the laboratory and the instruments

Six Sigma tools allow the laboratory to

- Identify the **RIGHT** method
- Select the **RIGHT** rules
- Run the **RIGHT** number of controls
- *(NEW! Run controls at the **RIGHT** frequency!)*
- Most importantly, enable the **RIGHT** patient outcomes



**ENTER YOUR
VIEWS IN OUR POC
QC SURVEY ON
JAMES.WESTGARD.
COM**



**QUESTIONS?
STEN@WESTGARD.COM**